

Cardiac Output in the Australian Bearded Lizard – Case Study



USCOM Case Study

Cardiac Output in the Australian Bearded Dragon

Rob Phillips¹, Julie Sloniger², Jennifer Hess³ Dr Martinez³

1. Department of Medicine, University of Queensland, Brisbane, Australia.
2. USCOM Ltd, Clinical Applications Specialist, Tucson, USA.
3. Texas A and M, Dept of Veterinary surgery, Austin, USA.

Background:

Circulation is an almost universal phenomenon in living organisms functioning to deliver the substrates for cell and tissue performance and removing the metabolic wastes from these biochemical reactions. Our understanding of circulation is substantially based on invasive assessments such as thermodilution, a method which presents a delayed, damped and time averaged representation of circulation, or single snapshot methods such as Fick. USCOM is a noninvasive ultrasound device specialized for beat to beat serial and accurate measurement of flow with very high sensitivity, and is currently used across multiple clinical and research applications in humans. The device is also being adopted in animal research and veterinary applications because of an increased focus on therapeutic veterinary haemodynamics and ethical animal experimentation. The device is currently being used to evaluate circulation in rats, dogs, pigs, sheep and foals. This study was to determine the feasibility of evaluating the central circulation of a reptile, the Black Bearded Lizard.

Method:

Two captive Bearded Dragons were examined using the USCOM device and central haemodynamics measured. The Lizards were held in the erect position (figure 1), and the ultrasound transducer placed approximately midline on the ventral aspect of the thorax and directed from a parasternal access toward the dorsum and the dominant flow profile was detected and optimized (figure 2). Velocity time flow profiles were



Uscom Limited
ABN 35 091 028 090
Level 7, 10 Loftus Street
Sydney NSW 2000 Australia
T +612 9247 4144 F +612 9247 8157
www.uscom.com.au

Cardiac Output in the Australian Bearded Lizard – Case Study



acquired and measured using the conventional ultrasound method while the USCOM algorithm was applied for calculation of flow volume. Flow profiles are demonstrated in figures 3 and 4. Cardiac Output (CO), stroke volume (SV), and heart rate (HR) were measured and values indexed to the mass of the reptiles for comparison. The lizards were returned to their pen after examination and were without obvious distress. Observations were made at 30.6°N of the Equator, at a temperature of 18°, on the 27th Julian day, at 15hrs past midnight.




Figure 1. Bearded Dragon undergoing USCOM examination during support in the erect position.



Cardiac Output in the Australian Bearded Lizard – Case Study



Figure 2. Bearded Dragon comfortable and at rest during a non-invasive USCOM examination via the parasternal acoustic access.

Results:

USCOM examinations




Figure 3. USCOM screen from lizard 1 (150 gms) demonstrating an SV of 0.77 cm³, a HR of 77 bpm and a CO of 50 cm³/min.

Cardiac Output in the Australian Bearded Lizard – Case Study



The measure of life.




Figure 4. USCOM screen from lizard 2 (250 gms) demonstrating an SV of 1.4 cm³, a HR of 67 bpm and a CO of 90 cm³/min.

Lizard 1 (150 gms): SV = 0.77 cm³, HR = 77 bpm, CO = 50 cm³/min.

Lizard 2 (250 gms): SV = 1.4 cm³, HR = 67 bpm, CO = 90 cm³/min.

CO indexed to mass = 333cm³/min/kg for lizard 1 and 360c m³/min/kg for lizard 2.

SV indexed to mass = 5.1 cm³/min/kg for lizard 1 and 5.6 cm³/min/kg for lizard 2

Discussion:

USCOM demonstrated easily accessible circulatory information in two small reptiles without the need for intervention or post examination euthanasia. The lizards showed no signs of distress throughout and resumed normal captive activities after the examination. While data on central Bearded Dragon haemodynamics is scant, the results measured by USCOM were not considered unreasonable and the mass indexed values of CO and SV showed close agreement for the two lizards. The flow volumes were derived from the USCOM neonatal algorithm and may not be transposable to the Lizards. However the values seem feasible, and even without the



Uscom Limited
ABN 35 091 028 090
Level 7, 10 Loftus Street
Sydney NSW 2000 Australia
T +612 9247 4144 F +612 9247 8157
www.uscom.com.au

Cardiac Output in the Australian Bearded Lizard – Case Study



algorithm the flow data will accurately demonstrate changes in flow. The impact of orthostasis was not tested in these examinations.

Conclusion:

It is a feasible using USCOM to non-invasively access central circulation in Australian Bearded Dragons. Validation of normal values remains problematic but the presented values do not seem exceptional and the mass indexed values remarkably consistent. USCOM may provide an ethical insight into the circulation of small reptiles.

Teaching Points:

USCOM provides insight into central circulation of many animals, including small reptiles.

USCOM is a noninvasive, simple and ethical method with multiple veterinary and zoologic applications including animal research.

References:

- * Knobloch K, Lichtenberg A, Winterhalter M, Rossner D, Pichlmaier M, Phillips R. Noninvasive Cardiac Output Determination by Two-Dimensional Independent Doppler During and After Cardiac Surgery. *Ann Thorac Surg* 2005;80:1479-84
- * Lester A Critchley, Zhi Y Peng, Benny S Fok, Jules Flach, Simon C Wong, Anna Lee, Robert A Phillips. Testing the reliability of a new Ultrasonic Cardiac Output Monitor, the USCOM, using aortic flow probes in anaesthetized dogs. *Anesth Analg* 2005;100:748-53
- * Phillips RA, Hood SG, Jacobson BM, Lichtenthal PR, Burstow DJ, West MJ, May CN. Measurement of CO by flow probe, USCOM and PAC in conscious sheep at rest and after dobutamine. 26th ISICEM 2006, *Critical Care* 2006; 10(Sup1): 138
- * Shelton G, Jones DR. The physiology of the alligator heart – The cardiac cycle. *Journal of Experimental Biology* 1991; 158: 539-564
- * Malvin GM, Hicks JW, Greene ER Central vascular flow patterns in the alligator *Alligator Mississippiensis*. *Am J of Physiol – regulatory, integrative and comparative physiology* 1995; 269 (5): R1133-R1139
- * Kik MJL, Mitchell MA. Reptile cardiology: A review of anatomy and physiology, diagnostic approaches, and clinical disease. *Seminars in Avian and Exotic Pet Medicine* 2005; 14 (1): 52-60
- * Frappell P, Schultz T, Christian K. Oxygen transfer during aerobic exercise in a varanid lizard *Varanus mertensi* is limited by the circulation. *Journal of Experimental Biology* 2002; 205 (17): 2725-2736

